

PTBP2 Rabbit mAb

Catalog No: #57155

Package Size: #57155-1 50ul #57155-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PTBP2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	PTBP2 Antibody detects endogenous levels of total PTBP2
Immunogen Description	A synthesized peptide derived from human PTBP2
Other Names	MIBP; nPTB; nPTB5; nPTB6; nPTB7; nPTB8; PTB; PTBLP; Ptbp2;
Accession No.	Q9UKA9
Calculated MW	55kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of PTBP2 expression in Neuro2a cell lysate.

Background

RNA-binding protein which binds to intronic polypyrimidine tracts and mediates negative regulation of exons splicing.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.