

ORM1 Conjugated Antibody

Catalog No: #C56212

Package Size: #C56212-AF350 100ul #C56212-AF405 100ul #C56212-AF488 100ul #C56212-AF555 100ul #C56212-AF594 100ul #C56212-AF647 100ul #C56212-AF680 100ul #C56212-AF750 100ul #C56212-Biotin 100ul #C56212-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ORM1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	ORM1 Antibody detects endogenous levels of total ORM1
Immunogen Description	A synthesized peptide derived from human ORM1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AGP 1; AGP A; AGP1; Alpha 1 acid glycoprotein; alpha-1-AGP; HEL S 153w; ORM; Orosomucoid 1;
Accession No.	Uniprot:P02763
Calculated MW	47kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Appears to function in modulating the activity of the immune system during the acute-phase reaction.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.