

## GSDMD (N terminal) Rabbit pAb

Catalog No: #55651



Package Size: #55651-1 50ul #55651-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	GSDMD (N terminal) Rabbit pAb
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB, IF/ICC
Species Reactivity	Human, Mouse, Rat
Immunogen Description	Recombinant fusion protein of human GSDMD(NP_079012.3).
Other Names	DF5L;DFNA5L;FKSG10;GSDMDC1;GSDMD
Calculated MW	53kDa
Storage	Store at -20°C. Avoid freeze / thaw cycles.   Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.

## Application Details

WB, 1:500 - 1:2000

IF/ICC, 1:50 - 1:200

## Images



Immunofluorescence analysis of A-431 cells using GSDMD (Full Length+N terminal) Rabbit pAb (A20197) at dilution of 1:200 (40x lens). Blue: DAPI for nuclear staining.

## Background

Gasdermin D is a member of the gasdermin family. Members of this family appear to play a role in regulation of epithelial proliferation. Gasdermin D has been suggested to act as a tumor suppressor. Alternatively spliced transcript variants have been described.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.