

GABAA R $\alpha$ 4 Polyclonal Antibody

Catalog No: #48630



Package Size: #48630-1 50ul #48630-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	GABAA R $\alpha$ 4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	GABAA R $\alpha$ 4 Polyclonal Antibody detects endogenous levels of GABAA R $\alpha$ 4 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human GABRA4.
Other Names	GABRA4; Gamma-aminobutyric acid receptor subunit alpha-4; GABA(A) receptor subunit alpha-4
Calculated MW	60kD
Concentration	1mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1:500 - 1:2000.

ELISA: 1:40000

## Images



Western Blot analysis of various cells using GABAA R $\alpha$ 4 Polyclonal Antibody diluted at 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.