

# ATP5G Rabbit mAb

Catalog No: #57225



Package Size: #57225-1 50ul #57225-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                                          |
|-----------------------|------------------------------------------------------------------------------------------|
| Product Name          | ATP5G Rabbit mAb                                                                         |
| Host Species          | Rabbit                                                                                   |
| Clonality             | Monoclonal                                                                               |
| Isotype               | Rabbit IgG                                                                               |
| Applications          | WB IHC IP                                                                                |
| Species Reactivity    | Human,Rat                                                                                |
| Immunogen Description | A synthetic peptide of human ATP5G1/G2/G3                                                |
| Other Names           | ATP synthase lipid-binding protein; ATP synthase membrane subunit c locus 1              |
| Accession No.         | Swiss-Prot:P05496/Q06055/P48201 Gene ID:516/517/518                                      |
| Calculated MW         | 14 kDa                                                                                   |
| SDS-PAGE MW           | 14 kDa                                                                                   |
| Formulation           | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA    |
| Storage               | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |


## Application Details

IP:1/20WB: 1/500-1/1000IHC: 1/50-1/100

## Images



Western blot analysis of ATP5G1/G2/G3 in Jurkat, C6 lysates using ATP5G antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using ATP5G1/G2/G3 antibody.High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.