

# Mouse YTH domain family protein 1 (YTHDF1) ELISA Kit

Catalog No: #EK20612

Package Size: #EK20612-1 48T #EK20612-2 96T



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse YTH domain family protein 1 (YTHDF1) ELISA Kit
Brief Description	Mouse YTH domain family protein 1 (YTHDF1) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

## Application Details

Detect Range:31.2-2000 pg/mL
Sensitivity:11.6pg/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 &mu;L
Assay Time:1-4.5h
Detection wavelength:450 nm
Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.