

Arrestin- β -1 Polyclonal Antibody

Catalog No: #48857



Package Size: #48857-1 50ul #48857-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Arrestin- β -1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Arrestin- β -1 Polyclonal Antibody detects endogenous levels of Arrestin- β -1 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human ARRB1.
Other Names	ARRB1; ARR1; Beta-arrestin-1; Arrestin beta-1
Calculated MW	50kD
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C 1 year

Application Details


WB 1:500 - 1:2000.

IHC 1:100 - 1:300.

IF 1:200 - 1:1000.

ELISA: 1:10000.

Images



Western blot analysis of lysates from Jurkat and HUVEC cells, using ARRB1 Antibody. The lane on the right is blocked with the synthesized peptide.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.