

Human Ig alpha-1 chain C region (IGHA1) ELISA Kit

Catalog No: #EK20653



Package Size: #EK20653-1 48T #EK20653-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Ig alpha-1 chain C region (IGHA1) ELISA Kit
Brief Description	Human Ig alpha-1 chain C region (IGHA1) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:1.56-100µg/mL

Sensitivity:0.52µg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.