

Human Ig alpha-2 chain C region (IGHA2) ELISA Kit

Catalog No: #EK20654



Package Size: #EK20654-1 48T #EK20654-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Ig alpha-2 chain C region (IGHA2) ELISA Kit
Brief Description	Human Ig alpha-2 chain C region (IGHA2) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:7.81-500ng/mL

Sensitivity:3.2ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.