

ASCL1 Polyclonal Antibody

Catalog No: #49035



Package Size: #49035-1 50ul #49035-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ASCL1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	ASCL1 Polyclonal Antibody detects endogenous levels of ASCL1 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human ASCL1.
Target Name	ASCL1
Other Names	ASCL1;ASH1;BHLHA46;HASH1;Achaete-scute homolog 1;ASH-1;hASH1;Class A basic helix-loop-helix protein 46;bHLHa46
Calculated MW	25kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500 - 1:2000.
IHC 1:100 - 1:300.
ELISA: 1:40000..
IF 1:50-200

Images



Western blot analysis of lysate from Jurkat cells, using ASCL1 antibody.



Immunohistochemistry analysis of ASCL1 antibody in paraffin-embedded human brain tissue.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.