

Human Symmetric dimethylarginine (SMDA) ELISA Kit

Catalog No: #EK20658



Package Size: #EK20658-1 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Symmetric dimethylarginine (SMDA) ELISA Kit
Brief Description	Human Symmetric dimethylarginine (SMDA) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range: 3.12 ng/mL ~ 200 ng/mL

Sensitivity: 1.56 ng/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time: 1-4.5 h

Detection wavelength: 450 nm

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.