

Rat Actin, aortic smooth muscle (ACTA2) ELISA Kit

Catalog No: #EK20673



Package Size: #EK20673-1 48T #EK20673-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Actin, aortic smooth muscle (ACTA2) ELISA Kit
Brief Description	Rat Actin, aortic smooth muscle (ACTA2) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:6.25-400 ng/mL

Sensitivity:2.87ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.