

# GAPDH(Internal) Monoclonal Antibody

Catalog No: #27031



Package Size: #27031-1 50ul #27031-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	GAPDH(Internal) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu P Rt Ms
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 72aa-248aa of human GAPDH (Internal).
Target Name	GAPDH
Other Names	glyceraldehyde-3-phosphate dehydrogenase □aging-associated gene 9 protein □peptidyl-cysteine S-nitrosylase GAPDH □glyceraldehyde 3-phosphate dehydrogenase
Accession No.	UniProt:P04406 NCBI Gene ID:2597
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

## Application Details

Predicted MW: 35kd  
Western blotting: 1:1000

## Images



1:1000 dilution of this antibody detected GAPDH on 40ug of cell lysates

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.