

Syndecan 1 Antibody Biotin Conjugated

Catalog No: #C02324B

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Syndecan 1 Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB ELISA IHC-P IHC-F Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 275-310 310 derived from human CD138
Conjugates	Biotin
Target Name	CD138 Syndecan 1
Other Names	SDC; CD138; SYND1; syndecan; Syndecan-1; SDC1
Accession No.	Swiss-Prot#P18827NCBI Gene ID6382
Cell Localization	Cytoplasm, Membrane bound
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

WB=1:500-2000

ELISA=1:5000-10000

IHC-P=1:100-500

IHC-F=1:100-500

Flow-Cyt=1&mu;g/Test IF=1:100-500

Background

Cell surface proteoglycan that bears both heparan sulfate and chondroitin sulfate and that links the cytoskeleton to the interstitial matrix.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.