STAT1(Phospho-Tyr701) Antibody

Catalog No: #11044

Package Size: #11044-1 50ul #11044-2 100ul



Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com

		4.5
	escri	ntion
$\boldsymbol{\nu}$	COUL	puon


Product Name	STAT1(Phospho-Tyr701) Antibody	
Host Species	Rabbit	
Clonality	Polyclonal	
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates.	
	Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho	
	specific antibodies were removed by chromatogramphy using non-phosphopeptide.	
Applications	WB IHC IP	
Species Reactivity	Human;Mouse;Rat	
Specificity	The antibody detects endogenous level of STAT1 only when phosphorylated at tyrosine 701.	
Immunogen Type	Peptide-KLH	
Immunogen Description	Peptide sequence around phosphorylation site of tyrosine 701 (T-G-Y(p)-I-K) derived from Human STAT1.	
Conjugates	Unconjugated	
Target Name	STAT1	
Modification	Phospho	
Other Names	CANDF7, ISGF-3, STAT91	
Accession No.	Swiss-Prot: P42224NCBI Protein: NP _009330.1	
SDS-PAGE MW	84,91kd	
Concentration	1.0mg/ml	
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02%	
	sodium azide and 50% glycerol.	
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.	

Application Details


Predicted MW: 84,91kd WB: 1:500~1:1000

IHC: 1:50~1:100IP: 1:50-200


Images



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using STAT1(Phospho-Tyr701) Antibody #11044(left) or the same antibody preincubated with blocking peptide(right).



Western blot analysis of extracts from MEF cells untreated or treated with interferon-I (IFNI) using STAT1 (Phospho-Tyr701) Antibody #11044.



Western blot analysis of extracts from Hela cells, treated with IFNa or calf intestinal phosphatase (CIP), using STAT1 (Phospho-Tyr701) Antibody #11044.

Background

Signal transducer and activator of transcription that mediates signaling by interferons (IFNs). Following type I IFN (IFN-a and IFN-beta) binding to cell surface receptors, Jak kinases (TYK2 and JAK1) are activated, leading to tyrosine phosphorylation of STAT1 and STAT2. The phosphorylated STATs dimerize, associate with ISGF3G/IRF-9 to form a complex termed ISGF3 transcription factor, that enters the nucleus. ISGF3 binds to the IFN stimulated response element (ISRE) to activate the transcription of interferon stimulated genes, which drive the cell in an antiviral state. In response to type II IFN (IFN-gamma), STAT1 is tyrosine- and serine-phosphorylated. It then forms a homodimer termed IFN-gamma-activated factor (GAF), migrates into the nucleus and binds to the IFN gamma activated sequence (GAS) to drive the expression of the target genes, inducing a cellular antiviral state.

Heim M H, (1999) J Recept Signal Transduct Res. 19: 75-120.

Durbin J E, et al. (1996) Cell. 84: 443-450.

Meraz M A, et al. (1996) Cell. 84: 431-442.

Wakao H, et al. (1994) EMBO J. 13: 2182-2191.

Demoulin J, B. et al. (1999) J Biol Chem. 274: 25855-258

Ihle J N, et al. (1994) Trends Biochem Sci. 19: 222-227.

Published Papers

C.-C.E. Lan, C.-S. Wu, S.-M. Huang el at., High-glucose environment reduces human b-defensin-2 expression in human keratinocytes: implications for

poor diabetic wound healing., British Association of Dermatologists, 166(6):1221n— C1229(2012)

PMID:22283836

Hirokazu Hara, Yoko Nakamura, Masayuki Ninomiya el at., Inhibitory effects of chalcone glycosides isolated from?Brassica rapa?L.

ı ε1• hidabeniı η and their synthetic derivatives on LPS-induced NO production in microglia., Bioorganic & Medicinal Chemistry, 19(18):5559-5568(2011)

PMID:21856162

Hsin-Chien Chen, Hsin-I Ma, Huey-Kang Sytwu el at., Neural Stem Cells Secrete Factors That Promote Auditory Cell Proliferation Via a Leukemia Inhibitory Factor Signaling Pathway, Journal of Neuroscience Research, 88:3308n— C3318(2010)

PMID:20882565

Jelke J. Fros, Wen Jun Liu, Natalie A. Prow el at., Chikungunya virus nonstructural protein 2 inhibits type I/II interferon-stimulated JAK-STAT signaling., Journal of Virology, 84(20)10877-10887(2010)

PMID:20686047

el at., Danzhi Jiangtang Capsule ameliorates kidney injury via inhibition of the JAK-STAT signaling pathway and increased antioxidant capacity in STZ-induced diabetic nephropathy rats.In Biosci Trends.On 2019 Jan 22 by Sun M, Bu W et al..PMID: 30606979, (2018)

PMID:30606979

el at., Chronic Inflammation In Response to Injury: Retention of Myeloid Cells in Injured Tissue Is Driven by Myeloid Cell IntrInsic Factors. In J Invest Dermatol on 2019 Jan 28 by Torbica T, Wicks K, et al..PMID:30703358, , (2019)

PMID:30703358

el at., A thiazole-derived oridonIn analogue exhibits antitumor activity by directly and allosterically InhibitIng STAT3. In J Biol Chem on 2019 Nov 15;by Shen X, Zhao L,et al..PMID:31594861, , (2019)

PMID:31594861

el at., The C-terminal Domain of Chikungunya Virus nsP2 Independently Governs Viral RNA Replication, Cytopathicity, and Inhibition of Interferon Signaling.In J Virol on 2013 Sep by Jelke J Fros, Erika van der Maten, et al..PMID:23864632, , (2013)

PMID:23864632

el at., High-glucose Environment Reduces Human ϵ° Y-defensin-2 Expression in Human Keratinocytes: Implications for Poor Diabetic Wound Healing.In

Br J Dermatol on 2012 Jun by C-C E Lan, C-S Wu,et al..PMID:22283836, , (2012)

PMID:22283836

el at., Neural stem cells secrete factors that promote auditory cell proliferation via a leukemia inhibitory factor signaling pathway. In J Neurosci Res on 2010 Nov 15 by Hsin-Chien Chen, Hsin-I Ma, et al..PMID: 20882565, , (2010)

PMID:20882565

el at., Chikungunya virus nonstructural protein 2 inhibits type I/II interferon-stimulated JAK-STAT signaling. In J Virol on 2010 Oct by Jelke J Fros, Wen Jun Liu, et al..PMID:

20686047, , (2010)

PMID:20686047

el at., Inhibitory effects of chalcone glycosides isolated from Brassica rapa L. I. ζ_1^l idabenii. ?and their synthetic derivatives on LPS-induced NO production in microglia. In Bioorg Med Chem on 2011 Sep 15 by Hirokazu Hara, Yoko Nakamura, et al..PMID:21856162, , (2011)

PMID:21856162

el at., DJC suppresses Advanced Glycation End Products-Induced JAK-STAT signaling and ROS in mesangial cells.In Evid Based Complement Alternat Med on 2017 May 29 by Min Sun, Yan Li ,et al..PMID: 28630633, , (2017)

PMID:28630633

el at., Inhibitory Effects of Chalcone Glycosides Isolated From Brassica Rapa L. 'Hidabeni' and Their Synthetic Derivatives on LPS-induced NO Production in Microglia.In Bioorg Med Chem on 2011 Sep 15 by Hirokazu Hara, Yoko Nakamura, et al..PMID: 21856162, , (2017)

PMID:21856162

el at., Newly Synthesized ι ζ ¬idabeniι ?Chalcone Derivatives Potently Suppress LPS-Induced NO Production via Inhibition of STAT1, but Not NF-I-• B, JNK, and p38, Pathways in Microglia.In Biol Pharm Bull On 2014 by Hirokazu Hara , Ryoko Ikeda et al..PMID: 24882415, , (2014)

PMID:24882415

el at., Molecular analysis of curcumin-induced polarization of murine RAW264. 7 macrophages.In J Cardiovasc Pharmacol on 2014 Jun by Fangyuan Chen , Ning Guo et al..PMID:24709638, , (2014)

PMID:24709638

el at., Expression of interferon effector gene SART1 correlates with interferon treatment response against hepatitis B infection. In Mediators Inflamm on 2016 by Yong Li, Chuanlong Zhu et al..PMID: 28077916, , (2016)

PMID:28077916

el at., ycobacterium tuberculosis EspB protein suppresses interferon-η1¬-induced autophagy in murine macrophages.In J Microbiol Immunol Infect on 2016 Dec by DanDan Huang, Lang Bao et al..PMID:25641595 , , (2016)

PMID:25641595

el at., Costunolide improved dextran sulfate sodium-induced acute ulcerative colitis in mice through NF-kB, STAT1/3, and Akt signaling pathways. In Int Immunopharmacol on 2020 Jul by Fan Xie, Hai Zhang, et al..PMID:32413737, , (2020)

PMID:32413737

Ighodaro Igbe;XiaoFei Shen;Wei Jiao;Zhe Qiang;Teng Deng;Sheng Li;WanLi Liu;HanWei Liu;GuoLin Zhang;Fei Wang el at., Dietary quercetin potentiates the antiproliferative effect of interferon-α in hepatocellular carcinoma cells through activation of JAK/STAT pathway signaling by inhibition of SHP2 phosphatase, , (2017)

PMID:29371942

Rui Yang;Meng-Yan Deng;Lu-Kui Yang;Guan-Zhe Wang;Jun Ma;Qiang Wen;Na Gao;Hai-Ling Qiao el at., Identification of Cytochrome P450 2E1 as a Novel Target in Neuroinflammation and Development of Its Inhibitor Q11 as a Treatment Strategy., , (2025)

PMID:40122152

Note: This product is for in vitro research use only and is not intended for use in humans or animals.