

VASP(Phospho-Ser239) Antibody

Catalog No: #11158



Package Size: #11158-1 50ul #11158-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	VASP(Phospho-Ser239) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB IHC
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous level of VASP only when phosphorylated at serine 239.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around phosphorylation site of serine 238 (K-V-S(p)-K-Q) derived from Human VASP.
Conjugates	Unconjugated
Target Name	VASP
Modification	Phospho
Accession No.	Swiss-Prot: P50552NCBI Protein: NP_003361.1
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 50kd


Western blotting: 1:500~1:1000

Immunohistochemistry: 1:50~1:100

Images



Western blot analysis of extracts from serum-treated 293 and PMA-treated C6 cells using VASP(Phospho-Ser239) Antibody #11158.



Immunohistochemical analysis of paraffin-embedded human tonsil carcinoma tissue using VASP(Phospho-Ser239) Antibody #11158(left) or the same antibody preincubated with blocking peptide(right).

Background

Ena/VASP proteins are actin-associated proteins involved in a range of processes dependent on cytoskeleton remodeling and cell polarity such as axon guidance, lamellipodial and filopodial dynamics, platelet activation and cell migration. VASP promotes actin filament elongation. It protects the barbed end of growing actin filaments against capping and increases the rate of actin polymerization in the presence of capping protein. VASP stimulates actin filament elongation by promoting the transfer of profilin-bound actin monomers onto the barbed end of growing actin filaments. Plays a role in actin-based mobility of *Listeria monocytogenes* in host cells. Regulates actin dynamics in platelets and plays an important role in regulating platelet aggregation.

Wang HG, et al.

Published Papers

el at., Influence of short-term refrigeration on collagen-dependent signalling mechanisms in stored platelets. In *Cell Signal* on 2024 Oct by Julia Zeller-Hahn, Marius Bittl,et al..PMID:39048036, , (2024)

[PMID:39048036](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.