

c-kit(phospho-Tyr936) Antibody

Catalog No: #11539

Package Size: #11539-1 50ul #11539-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	c-kit(phospho-Tyr936) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	IHC WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of c-kit only when phosphorylated at tyrosine 936.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around phosphorylation site of tyrosine 936 (H-I-Y(p)-S-N) derived from Human c-Kit.
Target Name	c-kit
Modification	Phospho
Other Names	KIT; CD117; SCFR; PBT;
Accession No.	Swiss-Prot: P10721NCBI Protein: NP_000213.1
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 120 145kd


Western blotting: 1:500~1:1000

Immunohistochemistry: 1:50~1:100

Images



Western blot analysis of extracts from C6 cell untreated or treated with serum using c-kit(phospho-Tyr936) Antibody #11539.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using c-kit(Phospho-Tyr936) Antibody #11539(left) or the same antibody preincubated with blocking peptide(right).

Background

This is the receptor for stem cell factor (mast cell growth factor). It has a tyrosine-protein kinase activity. Binding of the ligands leads to the autophosphorylation of KIT and its association with substrates such as phosphatidylinositol 3-kinase (PI3K)

Martin, F.H. et al. (1990) Cell 63, 203-11.

Blume-Jensen, P. et al. (2000) Nat Genet 24, 157-62.

Gommerman, J.L. et al. (1997) J Biol Chem 272, 30519-25.

Nocka, K. et al. (1990) EMBO J 9, 1805-13.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.