

## TBC1D4 (Phospho-Thr642) Polyclonal Antibody

Catalog No: #12385



Package Size: #12385-1 50ul #12385-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TBC1D4 (Phospho-Thr642) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Human;Mouse
Specificity	Phospho-TBC1D4 (T642) Polyclonal Antibody detects endogenous levels of TBC1D4 protein only when phosphorylated at T642.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human TBC1D4 around the phosphorylation site of T642.
Conjugates	Unconjugated
Target Name	TBC1D4
Modification	Phospho
Other Names	TBC1D4; AS160; KIAA0603; TBC1 domain family member 4; Akt substrate of 160 kDa; AS160
Accession No.	Swiss-Prot: O60343NCBI Gene ID: 9882
Target Species	human
SDS-PAGE MW	150kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1:500-2000

Immunohistochemistry: 1/100 - 1/300

Immunofluorescence: 1/200 - 1/1000

ELISA: 1/5000

## Images



Western blot analysis of L929 using p-TBC1D4 (T642) antibody.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.