

E1A-associated Protein p300 (Phospho-Ser1834) Antibody

Catalog No: #12453

Package Size: #12453 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E1A-associated Protein p300 (Phospho-Ser1834) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	E1A-associated Protein p300 (Phospho-Ser1834) Antibody detects endogenous levels of E1A-associated Protein p300 only when phosphorylated at Ser1834
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human E1A-associated Protein p300 (Phospho-Ser1834)
Conjugates	Unconjugated
Target Name	E1A-associated Protein p300
Modification	Phospho
Other Names	EP300, E1A binding protein p300, Histone acetyltransferase p300, KAT3B, p300 HAT, E1A-associated protein p300, E1A-binding protein, 300kD, p300, RSTS2
Accession No.	Swiss-Prot#: Q09472NCBI Gene ID: 2033
Target Species	human
Calculated MW	264kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Published Papers

el at., Pyruvate Kinase M2 Promotes the Activation of Dendritic Cells by Enhancing IL-12p35 Expression. In Cell Rep on 2020 May 26 by Xin Jin, Wenxia Zhang, et al..PMID:32460017, , (2020)

[PMID:32460017](https://pubmed.ncbi.nlm.nih.gov/32460017/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.