

JAK3 (Phospho-Tyr904) Antibody

Catalog No: #12511



Package Size: #12511-1 50ul #12511-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	JAK3 (Phospho-Tyr904) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	JAK3 (Phospho-Tyr904) Antibody detects endogenous levels of JAK3 only when phosphorylated at Tyr904
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human JAK3 (Phospho-Tyr904)
Conjugates	Unconjugated
Target Name	JAK3
Modification	Phospho
Other Names	JAK3 (Phospho-Tyr904) Antibody
Accession No.	Swiss-Prot#: P52333NCBI Gene ID: 3718
Target Species	human
Calculated MW	125kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images

Western Blot analysis of lysates of HepG2 UV treated and HepG2, using primary antibody at 1:1000 dilution.



Note: This product is for in vitro research use only and is not intended for use in humans or animals.