

IRF-7 (Phospho-Ser477) Antibody

Catalog No: #12763

Package Size: #12763 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IRF-7 (Phospho-Ser477) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	IRF-7 (Phospho-Ser477) Antibody detects endogenous levels of IRF-7 only when phosphorylated at Ser477
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human IRF-7 (Phospho-Ser477)
Conjugates	Unconjugated
Target Name	IRF-7
Modification	Phospho
Other Names	IRF7, Interferon regulatory factor 7, IRF-7, IRF7A, IRF-7H
Accession No.	Swiss-Prot#: Q92985NCBI Gene ID: 3665
Target Species	human
Calculated MW	54kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Published Papers

el at., Corrigendum: ScopoletIn Suppresses Activation of Dendritic Cells and Pathogenesis of Experimental Autoimmune Encephalomyelitis by Inhibiting NF-κB Signaling. In Front Pharmacol on 2019 Sep 13 by Zhang F, Zhang Y, et al.. PMID:31543821, , (2019)

[PMID:31543821](https://pubmed.ncbi.nlm.nih.gov/31543821/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.