

EZH2 (Phospho-Thr487) Antibody

Catalog No: #12820

Package Size: #12820-1 50ul #12820-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EZH2 (Phospho-Thr487) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purification using immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	EZH2 (phospho T487) Antibody detects endogenous levels of EZH2 only when phosphorylated at Thr487
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human EZH2 around the phosphorylation site of Thr487.
Target Name	EZH2 (Phospho-Thr487)
Modification	Phospho
Other Names	EZH2 (phospho T487) Antibody
Accession No.	Swiss-Prot#: Q15910 NCBI Gene ID: 2146
Target Species	human
Calculated MW	85 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Published Papers

el at., EZH2 K63-polyubiquitination affecting migration in extranodal natural killer/T-cell lymphomaIn Clin EpigeneticsOn2023 Nov 29byBoheng Li , Qidi Zhou et al..PMID:38031139 , (2023)

[PMID:38031139](https://pubmed.ncbi.nlm.nih.gov/38031139/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.