

GAP43 Antibody

Catalog No: #21273



Package Size: #21273-1 50ul #21273-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GAP43 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	WB IHC IF
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous level of total GAP43 protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa.38~43 (Q-A-S-F-R) derived from Human GAP43.
Conjugates	Unconjugated
Target Name	GAP43
Other Names	B-50; BASP2; NEUM; PP46; axonal membrane protein GAP-43
Accession No.	Swiss-Prot: P17677NCBI Protein: NP_001123536.1
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details


Predicted MW: 43kd

Western blotting: 1:500~1:1000


Immunohistochemistry: 1:50~1:100

Immunofluorescence: 1:100~1:200


Images



Western blot analysis of extracts from mouse brain tissue using GAP43(Ab-41) Antibody #21273 and the same antibody preincubated with blocking peptide.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using GAP43(Ab-41) Antibody #21273(left) or the same antibody preincubated with blocking peptide(right).



Immunofluorescence staining of methanol-fixed HeLa cells using GAP43(Ab-41) Antibody #21273.

Background

GAP43 encoded by this gene has been termed a 'growth' or 'plasticity' protein because it is expressed at high levels in neuronal growth cones during development and axonal regeneration. This protein is considered a crucial component of an effective regenerative response in the nervous system. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.

Rachael L. Neve, et.al. (1998) J. Neurosci; 18: 7757.

Yiping Shen, et.al. (2002) J. Neurosci; 22: 239.

Chantal Gamby, et.al. (1996) J. Biol. Chem; 271: 26698.

Published Papers

el at., Responses of PKCε to cardiac overloads on myocardial sympathetic innervation and NET expression.In Auton Neurosci. On 2018 Mar by Hu B, Zhang J et al..PMID:29195789, , (2018)

[PMID:29195789](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.