

BIM Antibody

Catalog No: #21280

Package Size: #21280-1 50ul #21280-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BIM Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	WB IHC
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous level of total BIM protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa.67~71 (P-A-S-P-G) derived from Human BIM.
Conjugates	Unconjugated
Target Name	BIM
Other Names	BCL2-like protein 11; BCL2L11; BOD; Bcl-2 related ovarian death protein;
Accession No.	Swiss-Prot: O43521NCBI Protein: NP_006529.1
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 23kd


Western blotting: 1:500~1:1000

Immunohistochemistry: 1:50~1:100

Images



Western blot analysis of extracts from HeLa cells using BIM(Ab-69) Antibody #21280.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using BIM(Ab-69) Antibody #21280(left) or the same antibody preincubated with blocking peptide(right).

Background

Induces apoptosis. Isoform BimL is more potent than isoform BimEL. Isoform Bim-a1, isoform Bim-a2 and isoform Bim-a3 induce apoptosis, although less potent than the isoforms BimEL, BimL and BimS. Isoform Bim-gamma induces apoptosis.

Fukazawa H, et al. (2004) Mol Cancer Ther ; 3(10): 1281-1288

Luciano F, et al. (2003) Oncogene ; 22(43): 6785-6793

Published Papers

el at., Elucidating the role of the FoxO3a transcription factor in the IGF-1-induced migration and invasion of uveal melanoma cancer cells. In Biomed Pharmacother on 2016 Dec by Fengxia Yan , Rifang Liao et al..PMID: 27881235, , (2016)

[PMID:27881235](#)

el at., The antisense lncRNA of TAB2 that prevents oxidative stress to enhance the follicular growth in mammals. In Commun Biol on 2024 Oct 2 by Nian Li, Bing Yun,et al..PMID:39358475, , (2024)

[PMID:39358475](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.