

LC3 antibody

Catalog No: #21400

Package Size: #21400-1 50ul #21400-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	LC3 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous level of total LC3 protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa. 112-116 (V-Y-A-S-Q) derived from Human LC3.
Conjugates	Unconjugated
Target Name	LC3
Other Names	LC3A; ATG8E; MAP1ALC3; MAP1BLC3;
Accession No.	Swiss-Prot: Q9H492NCBI Protein: NP_115903.1
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 19kd

Western blotting: 1:500~1:1000

Images



Western blot analysis of extract from K-562, HeLa, JK and 293 cells using LC3 Antibody #21400

Background

Probably involved in formation of autophagosomal vacuoles (autophagosomes).

Published Papers

el at., PIWI-Interacting RNA-004800 Is Regulated by S1P Receptor Signaling Pathway to Keep Myeloma Cell Survival.In Front Oncol on 2020 Apr 15
by Ma H, Wang H, et al..PMID: 32351883, , (2020)

[PMID:32351883](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.