

RNA silencing suppressor Antibody

Catalog No: #21411

Package Size: #21411-1 50ul #21411-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RNA silencing suppressor Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	WB
Species Reactivity	Virus
Specificity	The antibody detects endogenous level of total RNA silencing suppressor protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa.9~13 (D-A-R-E-Q) derived from Cucumber necrosis virus RNA silencing suppressor.
Target Name	RNA silencing suppressor
Accession No.	Swiss-Prot: P15184 NCBI Protein: NP_040957.1
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 20kd

Western blotting: 1:500~1:1000

Background

plays a role in initiation of gene silencing at the early stage of viral infection.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.