

RE57611p Antibody

Catalog No: #21416

Package Size: #21416-1 50ul #21416-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RE57611p Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	WB
Species Reactivity	Dm
Specificity	The antibody detects endogenous level of total RE57611 protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa.178~182 (K-R-Q-T-L) derived from Fruit fly RE57611p.
Target Name	RE57611p
Accession No.	Swiss-Prot: D3DMZ6 NCBI Protein: ADC54212
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 23kd

Western blotting: 1:500~1:1000

Background

CG8889 gene product from transcript CG8889-RA

Note: This product is for in vitro research use only and is not intended for use in humans or animals.