

Rex1(ZFP42) Antibody

Catalog No: #21661

Package Size: #21661-1 50ul #21661-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Rex1(ZFP42) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Rex1(ZFP42) protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa.10~14(K-T-R-H-Q) derived from Human Rex1(ZFP42).
Target Name	Rex1(ZFP42)
Other Names	REX1; ZNF754;
Accession No.	Swiss-Prot: Q96MM3NCBI Protein: NP_777560.2
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 35kd

Western blotting: 1:500~1:1000

Images



Western blot analysis of extracts from JK, HL-60 and Raji cells using Rex1(ZFP42) Antibody #21661.

Background

Involved in self-renewal property of ES cells. May be involved in transcriptional regulation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.