

Puma Antibody

Catalog No: #21688


Package Size: #21688-1 50ul #21688-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Puma Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide-KLH
Target Name	Puma
Other Names	BBC3, JFY1, JFY-1,
Concentration	1.0mg/ml

Images



Western blot analysis of extracts from 293 and MCF-7 cells using Puma Antibody #21688.

Published Papers

el at., Small molecule MIRA-1 induces in vitro and in vivo anti-myeloma activity and synergizes with current anti-myeloma agents. On Br J Cancer. On 2014 Apr 29 by Saha MN, Chen Y et al.. PMID:24691427, , (2014)

[PMID:24691427](https://pubmed.ncbi.nlm.nih.gov/24691427/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.