

DLST antibody

Catalog No: #22024

Package Size: #22024 100ul



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	DLST antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 188 and 423 of human DLST
Target Name	DLST
Accession No.	NCBI Gene ID: 1743NCBI mRNA#: NM_001933NCBI Protein#: NP_001924
Concentration	0.7mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 49kd
Western blotting: 1:500-1:3000
Immunohistochemistry: 1:50-1:500
Immunofluorescence: 1:100-1:200


Images




Sample (30 ug of whole cell lysate)
A: A431
10% SDS PAGE
Primary antibody diluted at 1: 1000



A:Immunohistochemical analysis of paraffin-embedded RT2 xenograft, using DLST antibody at 1: 500 dilution.
B:Immunohistochemical analysis of paraffin-embedded Cal27 xenograft, using DLST antibody at 1: 500 dilution.



Costained with Hoechst 33342



Immunofluorescence analysis of methanol-fixed A431, using DLST antibody at 1: 500 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.