

Ladinin 1 antibody

Catalog No: #22037

Package Size: #22037 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ladinin 1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 300 and 499 of Ladinin 1
Target Name	Ladinin 1
Accession No.	NCBI Gene ID: 3898NCBI mRNA#: BC009742NCBI Protein#: AAH09742
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 57kd

Immunohistochemistry: 1:100-1:250

Immunofluorescence: 1:100-1:200

Images




Sample (30 ug of whole cell lysate)


A: Raji

7.5% SDS PAGE


Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded Colon ca, using Ladinin 1 antibody at 1: 250 dilution.



Costained with Hoechst 33342



Immunofluorescence analysis of methanol-fixed A431, using Ladinin 1 antibody at 1: 500 dilution.

Background

The protein encoded by this gene may be an anchoring filament that is a component of basement membranes. It may contribute to the stability of the association of the epithelial layers with the underlying mesenchyme. [provided by RefSeq]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.