

Scinderin antibody

Catalog No: #22158

Package Size: #22158 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Scinderin antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 281 and 555 of human SCIN
Target Name	Scinderin
Accession No.	NCBI Gene ID: 85477NCBI mRNA#: NM_001112706NCBI Protein#: NP_001106177
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 81kd

Western blotting: 1:500-1:3000

Immunofluorescence: 1:100-1:200

Images




Sample (30 ug of whole cell lysate)

A: Hep G2

7.5% SDS PAGE

Primary antibody diluted at 1: 1000



Confocal immunofluorescence analysis (Olympus FV10i) of methanol-fixed HeLa, using SCIN antibody (Green) at 1: 500 dilution and alpha-tubulin antibody (Red) at 1: 2000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.