

## SHC4 antibody

Catalog No: #22280

Package Size: #22280 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	SHC4 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 352 and 609 of Human SHC4
Target Name	SHC4
Accession No.	NCBI Gene ID: 399694NCBI mRNA#: NM_203349NCBI Protein#: NP_976224
Concentration	0.5mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

Predicted MW: 69kd

Western blotting: 1:500-1:3000

## Images



Sample (30 ug of whole cell lysate)  
A: Molt-4  
7.5% SDS PAGE  
Primary antibody diluted at 1: 1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.