

PE anti-human CD140b

Catalog No: #248491

Package Size: #248491 25 Tests



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PE anti-human CD140b
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4A18A2
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Conjugates	Unconjugated
Other Names	Platelet-derived growth factor receptor, beta polypeptide, PDGFR1, PDGFRβ, PDGFRb, PDGF receptor bet
Formulation	0.2%BSAε 0.1%Proclin 300η£·ι 'Eη • ηO ε ²ζ¹AζA²oO PH7.2-7.4
Storage	4ε¹ oO ι• Ωε δΩ• ε-

Application Details

FC 5μl/ζ~¹

Note: This product is for in vitro research use only and is not intended for use in humans or animals.