

NANOG Antibody

Catalog No: #25045

Package Size: #25045 100ul




Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NANOG Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA WB IHC
Species Reactivity	Human;Mouse;Rat
Immunogen Type	Peptide
Immunogen Description	Raised against a 19 amino acid peptide near the center of human NANOG.
Conjugates	Unconjugated
Target Name	NANOG
Other Names	NANOG homeobox
Accession No.	EAW88651
Concentration	1mg/ml
Formulation	Supplied in PBS containing 0.02% sodium azide.
Storage	Can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Images



Western blot analysis of NANOG in human spleen tissue lysate with NANOG antibody at (A) 1 and (B) 2 ug/mL.



Immunohistochemistry of NANOG in human spleen tissue with NANOG antibody at 5 ug/mL.

Background

Expression of NANOG is required for the maintenance of pluripotency in epiblast and embryonic stem (ES) cells as well as for the ability to maintain ES self-renewal independently of LIF/Stat3. The role of NANOG in embryonic development suggested that it might be useful in the creation of stem cells that might be useful in cell replacement therapies in the treatment of several degenerative diseases. Artificial stem cells, termed induced pluripotent stem (iPS) cells, can be created by expressing POU5F1 (also known as Oct-4), another germline-specific transcription factor, and the transcription factors Sox2, Klf4 and Lin28 along with c-Myc in mouse fibroblasts. More recently, experiments have demonstrated that iPS cells could be generated using expression plasmids expressing NANOG, Sox2, Klf4 and c-Myc, eliminating the need for virus introduction, thereby addressing a safety concern for potential use of iPS cells in regenerative medicine.

Published Papers

et al., Wnt10B Is Critical for the Progression of Gastric Cancer. In Oncol Lett on 2017 Jun by Xiao-Dan Wu, Qing-Li Bie, et al.. PMID: 28599424, , (2017)
[PMID:28599424](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.