

RILP Antibody

Catalog No: #25472

Package Size: #25472 100ul




Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	RILP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA WB
Species Reactivity	Human;Mouse
Immunogen Type	Peptide
Immunogen Description	RILP antibody was raised against an 18 amino acid peptide near the center of human RILP.
Conjugates	Unconjugated
Target Name	RILP
Other Names	Rab-interacting lysosomal protein, PP10141
Accession No.	NP_113618
Concentration	1mg/ml
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

Images



Western blot analysis of RILP in HeLa cell lysate with RILP antibody at 1 ug/mL in (A) the absence and (B) the presence of blocking peptide.

Published Papers

Hai Qiu;Manting Zhang;Minchun Li;Chuanxiang Chen;Huijun Wang;Xia Yue et al., Methamphetamine and Methamphetamine-Induced Neuronal Exosomes Modulate the Activity of Rab7a via PTEN to Exert an Influence on the Disordered Autophagic Flux Induced in Neurons., (2025)
[PMID:40141286](https://pubmed.ncbi.nlm.nih.gov/40141286/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.