

Seasonal H1N1 Neuraminidase Monoclonal Antibody

Catalog No: #26046



Package Size: #26046 100ul


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Seasonal H1N1 Neuraminidase Monoclonal Antibody |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Clone No. | mAb (Clone 6F4C8) |
| Purification | Immunoaffinity chromatography purified IgG |
| Applications | ELISA |
| Species Reactivity | Virus |
| Immunogen Type | Peptide |
| Immunogen Description | Raised against a synthetic peptide from the seasonal influenza (H1N1) neuraminidase protein . |
| Conjugates | Unconjugated |
| Target Name | Seasonal H1N1 Neuraminidase |
| Other Names | Seasonal Influenza A H1N1 Neuraminidase (6F4C8): Influenza N1, flu N1, NA |
| Accession No. | ACA33620 |
| Concentration | 1mg/ml |
| Formulation | Supplied in PBS containing 0.02% sodium azide. |
| Storage | Can be stored at -20°C, stable for one year. |

Images



Seasonal influenza A N1 antibody (Cat. No. PM-5921) specifically recognizes seasonal (H1N1) N1, and does not cross-react with peptide corresponding to swine-origin influenza A (S-OIV, H1N1) N1 peptide, in ELISA.

Background

Influenza A virus has one of sixteen possible hemagglutinin (HA) surface proteins and one of nine possible neuraminidase (NA) surface proteins. In early 2009, a novel H1N1 swine-origin influenza (S-OIV) A virus was identified in specimens obtained from patients in Mexico and the United States. The genetic make-up of this swine flu virus is unlike any other: it is an H1N1 strain that combines a triple assortment first identified in 1998 including human, swine, and avian influenza with two new pig H3N2 virus genes from Eurasia, themselves of recent human origin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.