

PerCP-Cy5.5 anti-mouse CD185

Catalog No: #26844



Package Size: #26844-1 25T #26844-2 100T

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PerCP-Cy5.5 anti-mouse CD185
Host Species	Rat
Clonality	Monoclonal
Isotype	Rat IgG2b, k
Applications	FC
Species Reactivity	Mouse
Conjugates	Unconjugated
Storage	4 °C/0 °C 1 • 2ε δΩ • ε-

Application Details

FC 5μl/ζ~'

Note: This product is for in vitro research use only and is not intended for use in humans or animals.