

PE anti-mouse F4/80

Catalog No: #26848



Package Size: #26848-1 25T #26848-2 100T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	PE anti-mouse F4/80
Isotype	rat IgG2b
Applications	FC
Species Reactivity	Mouse
Conjugates	PE
Other Names	EMR1, Ly71
Formulation	0.2%BSAε 0.1%Proclin 300η£·ι Έη • ηΟ ε ²ζΙΑζΑ²οΟ PH7.2-7.4
Storage	4 ε ιοΟ ι• Ωε δΩ• ε-

Application Details

FC 5μ|ζι'

Note: This product is for in vitro research use only and is not intended for use in humans or animals.