

## JARID1C Monoclonal Antibody

Catalog No: #27002

Package Size: #27002-1 50ul #27002-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	JARID1C Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC
Species Reactivity	Human
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment(Human)
Conjugates	Unconjugated
Target Name	JARID1C
Other Names	lysine-specific demethylase 5C, □protein SmcX, □selected cDNA on X, □Smcy homolog, X-linked, □Smcx homolog, X chromosome, □histone demethylase JARID1C, □JmjC domain-containing protein SMCX, □Jumonji/ARID domain-containing protein 1C
Accession No.	UniProt:P41229 NCBI Gene ID:8242
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

## Application Details

Predicted MW: 220kd


Western blotting: 1:1000

## Images



1:1000 dilution of this antibody detected JARID1C on 40ug of cell lysates

1:150 dilution from a previous lot detected JARID1C in HeLa



---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.