

KIAA1967 Monoclonal Antibody

Catalog No: #27008

Package Size: #27008-1 50ul #27008-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	KIAA1967 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu Ms Rt P D
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 1aa-265aa of human KIAA1967(N-terminus).
Conjugates	Unconjugated
Target Name	KIAA1967
Other Names	p30 DBC protein DBC.1 <input type="checkbox"/> deleted in breast cancer 1 <input type="checkbox"/> deleted in breast cancer gene 1 protein
Accession No.	UniProt:Q8N163 NCBI Gene ID:57805
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details

Predicted MW: 120kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected KIAA1967 on 40ug of cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.