

## HDAC4 Monoclonal Antibody

Catalog No: #27009

Package Size: #27009-1 50ul #27009-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	HDAC4 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC
Species Reactivity	Human;Mouse;Rat
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 16aa-225aa of human HDAC4(N-terminus).
Conjugates	Unconjugated
Target Name	HDAC4
Other Names	histone deacetylase 4 □ histone deacetylase A
Accession No.	UniProt:P56524 NCBI Gene ID:9759
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

## Application Details

Predicted MW: 130kd


Western blotting: 1:1000

## Images



1:3000 dilution of this antibody detected HDAC4 on 40ug of cell lysates

1:150 dilution from a previous lot detected HDAC4 in Hela



---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.