

NPM1 Monoclonal Antibody

Catalog No: #27013



Package Size: #27013-1 50ul #27013-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	NPM1 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC
Species Reactivity	Hu Rb P Rt Ms
Immunogen Type	Recombinant protein
Immunogen Description	Full length recombinant protein(295aa)of human NPM1.
Conjugates	Unconjugated
Target Name	NPM1
Other Names	nucleophosmin □ numatrin □ nucleolar protein NO38 □ nucleophosmin/nucleoplasmin family, member
Accession No.	UniProt:P06748 NCBI Gene ID:4869
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details

Predicted MW: 38kd


Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected NPM1 on 40ug of cell lysates

1:200 dilution from a previous lot detected NPM1 in Hela



Note: This product is for in vitro research use only and is not intended for use in humans or animals.