

## Hsp70(N-term) Monoclonal Antibody

Catalog No: #27029

Package Size: #27029-1 50ul #27029-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	Hsp70(N-term) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu D P Ms Rt
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 96-289aa of human HSPA1A
Conjugates	Unconjugated
Target Name	Hsp70
Other Names	heat shock 70 kDa protein 1A/1B <input type="checkbox"/> HSP70-1/HSP70-2 <input type="checkbox"/> HSP70.1/HSP70.2 <input type="checkbox"/> heat shock 70kD protein 1A <input type="checkbox"/> heat shock-induced protein <input type="checkbox"/> heat shock 70 kDa protein 1/2 <input type="checkbox"/> dnaK-type molecular chaperone HSP70-1
Accession No.	UniProt:P08107 NCBI Gene ID:3303
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Storage	store at -20°C

## Application Details

Predicted MW: 70kd

Western blotting: 1:1000

## Images



1:1000 dilution of this antibody detected HSPA1A on 40ug of cell lysates

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.