

TRIM33 Monoclonal Antibody

Catalog No: #27032

Package Size: #27032-1 50ul #27032-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	TRIM33 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC
Species Reactivity	Hu D P Ms Rt
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 712aa-853aa of human TRIM33 (C-terminus).
Conjugates	Unconjugated
Target Name	TRIM33
Other Names	E3 ubiquitin-protein ligase TRIM33 □TIF1-gamma □protein Rfg7 □ectodermin homolog □RET-fused gene 7 protein □tripartite motif-containing 33 □transcriptional intermediary factor 1 gamma
Accession No.	UniProt:Q9UPN9 NCBI Gene ID:51592
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details


Predicted MW: 140kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected TRIM33(C-terminus) on 40ug of cell lysates



HeLa cells using anti- TRIM33 (C-terminus) antibody diluted 1:50

Note: This product is for in vitro research use only and is not intended for use in humans or animals.