

AKT1(C-term) Monoclonal Antibody

Catalog No: #27038



Package Size: #27038-1 50ul #27038-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	AKT1(C-term) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu Ms P Rt
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment of human AKT1 (C-terminus)
Conjugates	Unconjugated
Target Name	AKT1
Other Names	RAC-alpha serine/threonine-protein kinase <input type="checkbox"/> PKB alpha <input type="checkbox"/> RAC-PK-alpha <input type="checkbox"/> proto-oncogene c-Akt <input type="checkbox"/> protein kinase B alpha <input type="checkbox"/> rac protein kinase alpha
Accession No.	UniProt:P31749 NCBI Gene ID:207
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details

Predicted MW: 56kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected AKT1 on 40ug of cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.