

TFIIH(N-term) Monoclonal Antibody

Catalog No: #27053

Package Size: #27053-1 50ul #27053-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	TFIIH(N-term) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC
Species Reactivity	Hu Ms D P Rt
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 1aa-279aa of human GTF2H1 (N-terminus).
Conjugates	Unconjugated
Target Name	TFIIH
Other Names	general transcription factor IIH subunit 1 <input type="checkbox"/> BTF2 p62 <input type="checkbox"/> basic transcription factor 2 62 kDa subunit <input type="checkbox"/> general transcription factor IIH polypeptide 1 <input type="checkbox"/> TFIIH basal transcription factor complex p62 subunit
Accession No.	UniProt:P32780 NCBI Gene ID:2965
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Storage	store at -20°C

Application Details

Predicted MW: 62kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected GTF2H1 on 40ug of cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.