

## FAS Monoclonal Antibody

Catalog No: #27054



Package Size: #27054-1 50ul #27054-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	FAS Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Human
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 211aa-329aa of human CD95(C-terminus).
Conjugates	Unconjugated
Target Name	FAS
Other Names	tumor necrosis factor receptor superfamily member 6 □Fas AMA □FAS 827dupA □CD95 antigen □FASLG receptor □apoptosis antigen 1 □Delta Fas/APO-1/CD95 □APO-1 cell surface antigen □apoptosis-mediating surface antigen FAS □tumor nec
Accession No.	UniProt:P25445 NCBI Gene ID:355
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

## Application Details

Predicted MW: 45kd

Western blotting: 1:1000

## Images



1:1000 dilution of this antibody detected FAS on 40ug of cell lysates

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.