

Ikaros Mouse mAb

Catalog No: #27062

Package Size: #27062-1 50ul #27062-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Ikaros Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Purification	Purified by antigen-affinity chromatography.
Applications	IF/ICC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human IKZF1 expressed in E. Coli.
Conjugates	Unconjugated
Target Name	Ikaros
Other Names	DNA-binding protein Ikaros □CLL-associated antigen KW-6 □Ikaros (zinc finger protein) □ikaros family zinc finger protein 1 □lymphoid transcription factor LyF-1 □zinc finger protein, subfamily 1A, 1 (Ikaros
Accession No.	UniProt:Q13422NCBI Gene ID:10320
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	store at -20°C

Application Details


IF/ICC 1:50-1:200; FC 1:50-1:100

Images



Immunocytochemistry analysis of Ikaros in HeLa cells using Ikaros(C-term) Monoclonal Antibody (green), and DAPI(blue). Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin.

Flow Cytometry analysis of Jurkat cells using Ikaros(C-term)
Monoclonal Antibody (green) and negative control (red).



Note: This product is for in vitro research use only and is not intended for use in humans or animals.