

Rab2A(C-term) Monoclonal Antibody

Catalog No: #27064

Package Size: #27064-1 50ul #27064-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Rab2A(C-term) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC
Species Reactivity	Hu Ms D Rt P
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 39-206aa of human RAB2A(C-terminus).
Conjugates	Unconjugated
Target Name	Rab2A
Other Names	ras-related protein Rab-2A □ small GTP binding protein RAB2A □ RAB2, member RAS oncogene family
Accession No.	UniProt:P61019 NCBI Gene ID:5862
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Storage	store at -20°C

Application Details


Predicted MW: 24kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected RAB2A on 40ug of cell lysates



HeLa cells using anti- RAB2A(C-terminus)antibody diluted 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.